

Occhio

COFFEE-TRAK


Change from laser diffraction to sieve analysis with one click

Particle Characteristic Analyzer

7 microns to 5 mm

Dynamic Digital Image Analysis

The smart alternative to laser diffraction systems



**CLICK FOR
PRODUCT DEMONSTRATION**




COFFEE-TRAK

SIZE & SHAPE ANALYSIS FOR POWDERS & GRANULARS

Developed especilly for the finest grinds including Espresso & Turkish Coffees



SCHEDULE A **LIVE WEBINAR DEMONSTRATION**



COFFEE-TRAK: New particle size analyzer

OCCHIO SA: Since 2001, your partner in particle characterization


In 2001, an international and multidisciplinary **team of engineers**, created **OCCHIO SA**, in order to manufacture a **complete range of particle characteristic analyzers**, for particle applications beginning at **200 nm**. Whether it is for laboratory instrumentation, « AT LINE » or even « ON LINE » solutions, OCCHIO is prepared to be your partner in high-level powder characterization. OCCHIO brings you accuracy, profit and innovation.

Based on a combination of mechanical and vacuum dispersion, Coffee-Trak provides for fast and accurate, size and shape analysis. Because particle shapes can be identified and measured, conversion from sieving is seamless.
Tests take just seconds to complete



LEFT- About one teaspoon of ground coffee is fed in to the funnel.

RIGHT- Using vacuum dispersion, and a vibrating hopper, the sample flies laterally through the focal plane. Using monochromatic lighting, and special optics, precise images are captured for analysis.



Sieve analysis is labor intensive, messy, and subject to error, especially for non spherical particles.

Glass cleaning

Compressed air provides automatic cleaning of the glass surfaces.

No sticking particles during analysis

No need to clean optics after analysis




DRY & WET ANALYSIS

Full sieve distribution

Particle size and shape

Particle counting with Flowcell (WET option)

Kinetic mode



Applications: Sugar, Coffee, Chocolate, Fiber, Tobacco, Sand, Diamond, PTFE, Ceramics, Fertilizers...



COFFEE-TRAK: CALLISTO™

CALLISTO™ : Software for particle analysis


Since 2001, OCCHIO has developed revolutionary image analysis systems compatible with all particles types (powders, suspensions, emulsions and foams). These new systems, using the **powerful software CALLISTO™**, allow particle size measurements, shape analysis & particle counting. Now , more than 50 parameters are available ...

PARTICLE SIZE PARAMETERS

Area diameter

Inner diameter - Length

Width - Geodesic length - Mean diameter



PARTICLE SHAPE PARAMETERS

O. Aspect Ratio and Elongation: The Elongation is defined as $1-AspectRatio$ with $AspectRatio$ being the ratio between the width and the length of the particle

ISO Straightness: For very elongated particles, the straightness is ratio between the maximal distance and the Feret Length

ISO Compactness and ISO Roundness : Compactness and roundness are related to the degree to which the particle is similar to a disc. Roundness is less robust than compactness

ISO Circularity: Degree to which the particle (or its projection area) is similar to a circle

ISO Solidity: Solidity is the object area divided by the area enclosed by the convex hull (perfect to detect aggregates)

O. Bluntness : The Bluntness Index is the expression of a “maturity in the abrasion process”

O- Roughness: Amount of material to be removed from the shape before getting a smooth surface


Luminance: The luminance is the mean greyscale level of the particle

O-Porosity; Porosity estimator

SCATTERPLOT: Possibility of **having a representation** of the **particle size parameters** versus morphologic estimators. According to the particle size classes chosen by the operator, all plots are corresponding to particles. Operator can visualize each particle in clicking on each plot.

MICROSCOPE MODE: Display of each parameters and picture for selected particle on Scatterplot or live pictures. Operator can save picture for individual particle in Bitmap format.

(Particle size & shape distribution, trends, Excel export, sieves correlation...)



Coffee-Trak Options

Liquid Analysis

Autosampler


Do not hesitate to contact us for other idea...



| Model | COFFEE-TRAK |
|---|---|
| Particle size range | 7 microns – 5000 microns |
| Measurement time | 30 seconds to 10 minutes (sample depending) |
| Parameters | <p>ISO 9276-6; 7; 8 norms + OCCHIO parameters</p> <p>SIZE: ISO Area diameter; ISO Inner diameter; Mean diameter; Perimeter diameter; Crofton diameter; Half Crofton diameter; Width; Length; Ellipse Width; Ellipse Length; ISO Max Distance; ISO Geodesic Length;</p> <p>SHAPE: Occhio Bluntness; Occhio Roughness; Elongation; ISO Aspect Ratio; Ellipsoid Elongation; Ellipsoid Roundness; Ellipse Ratio; ISO Eccentricity</p> <p>ISO Straightness; ISO Roundness; ISO Compactness; ISO Extent; ISO Solidity; Convexity; ISO Circularity, Luminance mean</p> <p>Luminance var.</p> |
| Dimensions and weight | 89.5(L) x 47 (l) x (H)38cm, 28 kg |
| Disperser | Vibrating feeder + vacuum principle |
| Optics and imaging device | Camera 5Mpixels Gigabit Ethernet |
| Validation | IQ,OQ,PQ in option |
| Option | Wet module, autosampler |
| Computer specifications | Windows 7, Intel Core i5-650 @3.2GHz, 4MB cache ; 4 GB @ 1156MHz , HD 500GB Thanks to contact OCCHIO to check computer specifications |
| Power supply | 110-240 V 50/60 Hz |
| Working conditions | Temperature 5 °C – 45 °C, HR 35-80 % |
| Images format | Bitmap |
| Data storage | ‘.oph’ binary Occhio files format contains: Full size distribution values Shape and size percentiles Outline and greyscale levels of each particle |
| Statistics tools (distribution in number & volume) | Acquisition info (short overview of the used SOP) ; Size distribution; Size percentiles; Shape percentiles; Shape distribution; Mean shape by size; 2D scatter-plot (fully selectable particles map); 3D scatter-plot (include animation); Percentiles sample images; Sample images (BMP exportable format); Id card for each particle (BMP exportable format); Morphological and size filtering procedure ; Raw data export (text format); Table distribution export (text format); Table distribution and percentile export (Excel format); Automatic or custom reporting; Full image export (bmp format); Single particle image export (bmp format); Figure and graph export(bmp format) |



Reference code: OCC164 Coffee-Trak



Technical specifications


Working conditions

| | Description |
|---------------------|--|
| Working temperature | Temperature 5°C – 45°C Humidity 35% - 80% non condensing |
| Power Supply | 110 or 220 Vac 50-60Hz (auto-switching is not available please ask to Occhio to set power supply) Computer (supplied by Occhio) |

| | Description |
|------------------|---------------------------------------|
| Processor | Intel Core i7-2600 @3.4GHz, 4MB cache |
| Ram | 4 GB @ 1600MHz |
| Hard Disk | 500MB |
| Display | LCD, FullHD, 22" |
| Mouse, keyboard | USB (English) |
| Operating system | Windows Seven professional |

Optics and imaging device

| | Description |
|----------------------|--|
| Standard camera type | Camera 5Mpixels Gigabit Ethernet 2/3" interline progressive scan CCD |
| Camera resolution | 5.0 Millions of pixels 2448 x 2050 pixels |
| Pixel size | 3.45 µm |




| | |
|---|--|
| Lens type | Telecentric lens |
| Lens resolution (standard lens) (lens code 164-120-R1) | Calibration : 10 µm/Pixel |
| Field of view (standard lens) (lens code 164-120-R1) | 24480 x 20500µm @10µm/pixel |
| Lens resolution (Optional lens) (lens code 164-121-R1) | Calibration : 7 µm/Pixel |
| Field of view (Optional lens) (lens code 164-121-R1) | 17136 x 14350µm @7µm/pixel |
| Light source | Back light collimated monochromatic light source (wavelength 440nm) |
| Light wavelength | 440 nm (blue light) |

Dimensions and weight

| | Description |
|------------|---|
| Dimensions | 85 x 64 x 47 cm (Width x Height x Deep) |
| Weight | 52 kg |

Starting kit parts (these parts are included in the packing box at the delivery)

| Part number | Description | Quantity |
|---|---|----------|
| 164-081-R1  | Glass plate | 2 |
| 999-0014-R1* | Vacuum cleaning (220VAC-240VAC) | 1 |
| 999-0015-R1* | Vacuum cleaning (110VAC) | 1 |
| 999-0015-R1 | Vacuum cleaning filter | 1 |
| 164-500-R1 | Communication cable A | 1 |
| 164-501-R1 | Communication cable B | 1 |
| 999-0003-R1 | Power supply cable North America | 3 |
| 999-0004-R1 | Power supply cable Europe | 3 |
| 999-0013-R1** | Computer + LCD, FullHD, 21.5" + Mouse + Keyboard US | 1 |
| 999-0010-R1 | Keyboard USB (FR) | 1 |
| 999-1016-R1 | Calibrated quartz sand (Fr B 600µm 1180µm) Net weight 500g | 1 bottle |


Option 164108 external wet dispersion module

| | Description |
|--------------------------|---------------|
| Dispersion beaker volume | max 600 ml |
| Flow type | Circulation |
| Cell thickness | 1mm; 3mm; 6mm |

Coffee-Trak short instrument overview

Based upon a combination of mechanical and vacuum dispersion, the Coffee-Trak provides fast and accurate size and shape analysis of sieveable powders. The instrument combines a high quality imaging system with a robust mechanical design for at-line and on-line process control.

The powder is dispersed by means of a vacuum system, different layers of air flow inside the pipe are maintained in laminar conditions. The area of image capture is close to the powder entry point to avoid differential speed. The system uses a unique combination of telecentric lens and collimated monochromatic light; therefore all particles are within the focus region.



Sampling

The sample is fed directly onto the vibrating feeder. Using a standard operating procedure, the analysis duration and cleaning time are automatically set, and the measurement starts. Within a few seconds, results are displayed and a database is generated.

This instrument is able to analyze a large quantity of sample in just a few minutes due to the high-speed air flow.



All Occhio instruments are based upon specific optical technology, using high quality lenses with low distortion and mounted on a precision, robust mechanical system suitable for industrial working conditions.

Sample analysis

| | |
|---------------------------------------|---|
| Model | Occhio Coffee-Trak |
| Sample dispersion | Vibrating feeder combined with air flow |
| Sample particles size range | From 7 µm to 5 mm |
| Time measurement | 1-3 minutes (sample et quantity dependent) |
| Sample analysis | Size distribution cumulate and proportional curve Number distribution or volume weighted distribution Sieves correlation (according with OCCHIO) |
| Standard Operating Procedure includes | Maximum number of particles Control on vibrating feeder (initial speed, acceleration, max) Particles per picture Light intensity calibration Background calibration Creation of a particle database Image storage Filtering procedure Automatic report generation |

Software main features

| | |
|--|--|
| Model | Callisto Software |
| Size parameters (Iso 9276-6; 7; 8) All the size parameters are displayable or not according with the customer setting preference | ISO Area diameter ISO Inner diameter Mean diameter Perimeter diameter Crofton diameter Half Crofton diameter Width Length Ellipse Width Ellipse Length ISO Max Distance ISO Geodesic Length |
| Shape parameters (Iso 9276-6; 7; 8) All the shape parameters are displayable or not according with the customer setting preference | Occhio Bluntness Occhio Roughness Elongation ISO Aspect Ratio Ellipsoid Elongation Ellipsoid Roundness Ellipse Ratio ISO Eccentricity ISO Straightness ISO Roundness |



| | |
|---|--|
| | <p>ISO Compactness ISO Extent ISO Solidity Convexity ISO Circularity Luminance mean Luminance var. Porosity</p> |
| Advanced shape parameters | Developed in function of customer specifications |
| Image format | Bitmap |
| Data storage | <p>`.oph` binary Occhio files format contains: Full size distribution values Shape and size percentiles Outline and greyscale levels of each particle</p> |
| Data comparisons | Open and compare more analysis on the same plots include `trends graphic` |
| Plots and figure (By number or volume weighted values) | <p>Acquisition info (short overview of the used SOP) Size distribution Size percentiles Shape percentiles Shape distribution Mean shape by size 2D scatter-plot (fully selectable particles map) 3D scatter-plot (include animation) Percentiles sample images Sample images (BMP exportable format) Id card for each particle (BMP exportable format)</p> |
| Statistics tools | Morphological and size filtering procedure |
| Reporting and data export | <p>Raw data export (text format) Table distribution export (text format) Table distribution and percentile export (Excel format) Automatic or custom reporting Full image export (bmp format) Single particle image export (bmp format) Figure and graph export (bmp format)</p> |
| Microscope mode pane | Real time images acquisition without analyzing |
| External images analysis | Possibility to analyse images coming from others acquisition sources (bmp greyscale format) |

OCCHIO SA
4 rue des chasseurs ardennais BELGIUM Tel :+32 43729330 Fax : +32 43652346 info@occhio.be
www.occhio.be