

# **Sinar SP 6600 Moisture Analyser**

- Easy to use no sampling, weighing or grinding, simply insert and read
- Fast instant digital readout of moisture and temperature
- Robust rugged construction with no moving parts
- Accurate many readings can be taken in a few minutes to give a truly representative picture of moisture in bags or in bulk



#### **CLICK FOR QUOTATIONS & PRICING**

## **Product Description**

The SackPro allows the analysis of any commodity stored in bags without even having to take a sample, simply insert the probe to measure in seconds what traditional methods would take minutes to achieve. This makes the routine analysis of every sack a reality, giving you confidence that your product is exactly right without having to waste time taking samples back to the office or lab for analysis.

The SackPro comes calibrated for seven commodities which can be chosen from the large catalogue of calibrations that we have built up and refined over the past 20 years. Changing calibrations at a later date can also easily be achieved using the built in communication port. This means that however specialised the application the SackPro can deliver fast accurate results.

The SackPro has no moving parts and its mild steel polyester powder coated shaft is rugged enough to withstand even the toughest of operating conditions allowing results to be obtained wherever they are needed.

#### **Product Features**

Accuracy: typically  $\pm 0.5 \%$  mc and  $\pm 0.1 \%$ 

Measurement Range: 0 to 35 % mc dependant on application, -20 to 60  $^{\circ}$ C

Operating Environment: Dimensions:

Shaft 1 m, display box 20x15x3 cm

Weight: 1 kg

4 x 1.5V AA cells

Power Supply: Measurement principles:

Moisture - temperature compensated electric field Temperature - Platinum resistance detector

Printer output: RS232C, 300 or 4800 baud

Resolution: Moisture 0.1% Temperature 0.1℃

## **Supplier Address**

QAQC LAB / COFFEE LABORATORY 589 Rappahannock Drive White Stone Va 22578 United States

Tel (804) 435-5522

Toll Free US (866) 244-1578

Fax: (703) 991-7133

sales@coffeelabequipment.com